

```

! -----
! CYCLIC VOLTAMMETRY (CV)
! This is a test program for UWED to make basic CV recordings
! Written by Alar Ainla, 2016–2018. Whitesides Group, Harvard.
! -----

! -----
! OPERATION PARAMETERS
! -----

upperVoltage=0.5      ! [V] This is upper voltage value
lowerVoltage=-0.5     ! [V] This is lower voltage value
numberOfCycles=3      ! This is number of full cycles
scanRate=100          ! [mV/s] This is the scan rate used
timeStep=25           ! [ms] time step between measurements divided by 2.
! For example timesStep=25 would correspond to 50ms step between measurements

! -----
! PROGRAM CODE
! -----

! Here we lock working electrode potential W[d]=22500 W[V]=1.5031
! And only vary the reference electrode potential
! ADC[d]=(ADC[V])*2.01513E4+2.38491E2
! ADC[V]=(ADC[d]-2.38491E2)/2.01513E4
! Curr[uA]=0.00608*(dADC-30694)-0.04      ! Working electrode current
! (W-R) [V]=R[d]*(-6.64898E-5)+1.50310    ! Working electrode potential vs ref
Rslope=1.5040E4       ! This is number of digits per volt when Ref is changed
prp=0

! Calculate program data structure
DIM initW AS LONG
initW=22500           ! This is now constant and not calculated
DIM initR AS LONG
initR=22606           ! This should start as close to zero as possible
nrScans=2*numberOfCycles+1

DIM StepHeightAbs AS INTEGER
nrSteps1=1000000*(upperVoltage)/(scanRate*timeStep*2)
nrSteps2=1000000*(upperVoltage-lowerVoltage)/(scanRate*timeStep*2)
nrSteps3=nrSteps2-nrSteps1
stepHeightAbs=ABS((upperVoltage-lowerVoltage)*Rslope)/nrSteps2
DIM stepR(nrScans) AS INTEGER
DIM nrSteps(nrScans) AS INTEGER
DIM nrDataPoints AS LONG
nrDataPoints=2*numberOfCycles*(nrSteps2+2)
FOR i=1 TO numberOfCycles
    stepR(i*2)=stepHeightAbs
    stepR(i*2+1)=(-1)*stepHeightAbs
    nrSteps(i*2)=nrSteps2
    nrSteps(i*2+1)=nrSteps2
NEXT i
stepR(1)=(-1)*stepHeightAbs
nrSteps(1)=nrSteps1
nrSteps(nrScans)=nrSteps3
! Output these to console for feedback
PRINT "InitW"
PRINT initW

```

```

PRINT "InitR"
PRINT initR
PRINT "stepR"
PRINT stepR
PRINT "nrSteps"
PRINT nrSteps
PRINT "nrScans"
PRINT nrScans
PRINT "nrDataPoints"
PRINT nrDataPoints

! Resulting data storage
DIM Vw(nrDataPoints) AS DOUBLE ! Potential in volts
DIM IuA(nrDataPoint) AS DOUBLE ! Electrode current in uA
DIM dWe(nrDataPoints) AS LONG ! Working electrode as set
DIM dRe(nrDataPoints) AS LONG ! Reference electrode as set
DIM dADC(nrDataPoints) AS LONG ! ADC value
DIM Ts(nrDataPoints) AS LONG ! Time step
DIM Ni(nrDataPoints) AS INTEGER ! Index of point in scan
DIM DPi AS INTEGER ! Data point counter
DIM SCi AS INTEGER ! Scan counter
DIM ScanRunning AS INTEGER ! if 0 does not run, if 1 running
ScanRunning=0
DIM lastIndex AS INTEGER ! Last receiveing Index, to avoid double reading
lastIndex=0

! -- GUI --
! - Initialize GUI -
DIM startButton AS Button ! This button starts the measurement
DIM myMess AS Label ! This is text message area
DIM p1 AS Plot ! This is the result plot
DIM pCV AS PlotPoint ! This is one dataset
DIM myCV(nrDataPoints,2)
FOR i=1 TO nrDataPoints
    myCV(i,1)=0
    myCV(i,2)=0
NEXT i
p1=Graphics.newPlot
pCV=p1.newPlot(myCV)
p1.setView(lowerVoltage,-100,upperVoltage,100,0)
p1.setXaxisLabel("Potential Working vs Ref (V)")
p1.setYaxisLabel("Current working electrode (µA)")
p1.setTitle("Electrochemical CV with UWED")
startButton=Graphics.newButton(10,10)
startButton.setTitle("Start")
myMess=Graphics.newLabel(110,10,200,20)
myMess.setText("....")
System.ShowGraphics(1)

! -- BLE --
print "Hello BLE"
DIM bm AS BLEPeripheral
ble.startble
DIM uuid(0) AS STRING
DIM services(1) AS STRING
DIM bleWriteChar AS BLECharacteristic

ble.startscan(uuid)

```

END

```
! -----  
! FUNCTIONS  
! -----
```

SUB InitializeUWED()

DPi=0 ! No data points yet

SCi=0 ! Scan counter

! Make graph empty

FOR i=1 TO nrDataPoints

myCV(i,1)=0

myCV(i,2)=0

NEXT i

! Start setting up

sendToUWED("A("&STR(initR)&")")

sendToUWED("B("&STR(initW)&")")

sendToUWED("D(3)")

sendToUWED("C(1)")

END SUB

SUB InitializeScan()

SCi=SCi+1 ! Increase scan counter

IF SCi>nrScans THEN ! If scan counter is larger than number of scans, then DONE

sendToUWED("C(0)") ! Put UWED to potentiometric (off state)

! Save data

OPEN "CVData.txt" FOR OUTPUT AS #1

FOR i=1 TO DPi

PRINT #1, STR(dRe(i))&" "&STR(DWe(i))&" "&STR(dADC(i))&
" "&STR(Ts(i))&" "&STR(Ni(i))

NEXT i

CLOSE #1

! Send now data over email to server

DIM em AS Email

DIM em_mess AS STRING

em_mess="This is data from CV measurement from iPhone"

em=System.newEmail("email@emailserver.com","CV Data", em_mess)

em.addAttachment("CVData.txt", "text/plain")

em.send

myMess.setText("DONE")

ScanRunning=0

ELSE ! Initialize a new scan

! Increase same amount both steps

sendToUWED("G("&STR(stepR(SCi))&")")

sendToUWED("I("&STR(stepR(SCi))&")")

sendToUWED("L("&STR(nrSteps(SCi))&")")

sendToUWED("M()") ! Run the sequence

ScanRunning=1

END IF

END SUB

! - Process incoming data -

SUB procData(input_array() AS INTEGER)

PRINT "Index: "&CHR(input_array(1)) ! Feedback, so we can see continuity

! First check that there is no duplication

IF(input_array(1)=lastIndex) THEN

GOTO OutSub ! If same as last time do nothing

```

END IF
lastIndex=input_array(1) ! Otherwise put it same
! Now look if new data came M character on second position
IF(input_array(2)=ASC("M"))THEN ! New data came
    IF ScanRunning=1 THEN ! Only if scan is running there is something to come
        if DPi=20 and prp=0 then
            DPi=0
            prp=1
        end if
        DPi=DPi+1
        dRe(DPi)=ConvertToNumberBin(input_array,3,4)
        dADC(DPi)=ConvertToNumberBin(input_array,5,6)
        DPi=DPi+1
        dRe(DPi)=ConvertToNumberBin(input_array,7,8)
        dADC(DPi)=ConvertToNumberBin(input_array,9,10)
        Ts(DPi)=ConvertToNumberBin(input_array,11,14)
        Ts(DPi-1)=Ts(DPi)
        dWe(DPi)=ConvertToNumberBin(input_array,15,16)
        dWe(DPi-1)=dWe(DPi)
        Ni(DPi)=ConvertToNumberBin(input_array,17,18)
        Ni(DPi-1)=Ni(DPi)
        ! Do the math for visualization
        ! ODD
        VoltW=dRe(DPi-1)*(-6.64898E-5)+1.50310 ! Working electrode pot vs ref
        Curr=0.00608*(dADC(DPi-1)-30694)-0.04 ! Current in uA
        myCV(DPi-1,1)=VoltW
        myCV(DPi-1,2)=Curr
        ! EVEN
        VoltW=dRe(DPi)*(-6.64898E-5)+1.50310 ! Working electrode pot vs ref
        Curr=0.00608*(dADC(DPi)-30694)-0.04 ! Current in uA
        myCV(DPi,1)=VoltW
        myCV(DPi,2)=Curr
        ! Show graph
        pCV.setPoints(myCV)
        p1.setView(lowerVoltage,minY(myCV),upperVoltage,maxY(myCV),0)
    END IF
END IF
! If end of the scan start new
IF (input_array(2)=ASC("m")) THEN ! End of scan
    InitializeScan ! Start next scan
END IF
! In other cases do currently nothing
OutSub:
END SUB

```

! - Convert to Number from String source - currently not used

```

FUNCTION ConvertToNumberStr(input_array() AS INTEGER, _
    begi AS INTEGER, _
    endi AS INTEGER) AS LONG

```

```

    DIM nets AS STRING
    nets=""
    FOR i=begi TO endi
        nets=nets & chr(input_array(i))
    NEXT i
    ConvertToNumberStr=VAL(nets)
END FUNCTION

```

! - Convert to Number from Binary source - used in data receiving

```

FUNCTION ConvertToNumberBin(input_array() AS INTEGER, _
                           begi AS INTEGER, _
                           endi AS INTEGER) AS LONG

    DIM sum AS LONG
    sum=input_array(begi)
    FOR i=begi+1 TO endi
        sum=sum*256
        sum=sum+input_array(i)
    NEXT i
    ConvertToNumberBin=sum
END FUNCTION

! -----
! EVENTS
! -----

! When button is clicked
SUB touchUpInside(ctrl AS Button, when AS DOUBLE)
    IF ctrl=startButton THEN ! Send now
        myMess.setText("Start the scan")
        InitializeUWED
        InitializeScan
    END IF
END SUB

! BLE EVENTS and FUNCTIONS
! - Send command to UWED
SUB sendToUWED(inp AS STRING)
    DIM ax(LEN(inp)) AS INTEGER
    FOR i=1 TO LEN(INP)
        ax(i)=ASC(MID(inp,i,1))
    NEXT
    bm.writeCharacteristic(bleWriteChar,ax)
END SUB

! - First when device is discovered
SUB BLEDiscoveredPeripheral (    time AS DOUBLE, _
                               peripheral AS BLEPeripheral, _
                               services() AS STRING, _
                               advertisements(,) AS STRING, _
                               rssi)

    PRINT "Device Found";
    IF peripheral.bleName="UWED" THEN ! Right device found
        bm=peripheral
        ble.connect(bm)
        ble.stopScan
        print "Device is UWED"
    END IF
END SUB

! - When information about the device is fetched
SUB BLEPeripheralInfo(    time AS DOUBLE, _
                        peripheral AS BLEPeripheral, _
                        kind AS INTEGER, _
                        message AS STRING, _
                        err AS LONG)

    IF kind=1 THEN ! Connection completed

```

```

        peripheral.discoverServices(uuid)
ELSE IF kind=2 OR kind=3 THEN ! Connection lost
    ble.connect(bm)
ELSE IF kind=4 THEN ! Service dound
    DIM avServ(1) AS bleservice
    avServ=peripheral.services
    FOR a=1 TO UBOUND(services,1)
        FOR a=1 TO UBOUND(services,1)
            IF avServ(a).uuid="2220" THEN ! If right service id has been found
                peripheral.discoverCharacteristics(uuid, avServ(a))
                print "Service discovered"
            END IF
        NEXT
    NEXT
NEXT
END IF
END SUB

```

```

! - When information about service is obtained
SUB BLEServiceInfo( time AS DOUBLE, _
                    peripheral AS BLEPeripheral, _
                    service AS BLEService, _
                    kind AS INTEGER, _
                    message AS STRING, _
                    err AS LONG)

IF kind=1 THEN
    DIM chx(1) AS blecharacteristic
    chx=service.characteristics
    FOR i=1 TO UBOUND(chx,1)
        IF service.uuid="2220" THEN
            IF chx(i).uuid="2221" THEN
                peripheral.setNotify(chx(i),1) ! Notify if changes
                print "Read characteristic found"
            ELSE IF chx(i).uuid="2222" THEN
                bleWriteChar=chx(i)
                print "Write characteristic found"
            END IF
            myMess.setText("BLE CONNECTED!")
        END IF
    NEXT
END IF
END SUB

```

```

! - If BLE Characteristic info comes
SUB BLECharacteristicInfo( time AS DOUBLE, _
                           peripheral AS BLEPeripheral, _
                           characteristic AS BLECharacteristic, _
                           kind AS INTEGER, _
                           message AS STRING, _
                           err AS LONG)

IF kind=2 THEN
    IF characteristic.uuid="2221" THEN
        DIM value(1) AS INTEGER
        value=characteristic.value
        procData(value)
    END IF
END IF
END SUB

```

! -----
! HELPING FUNCTIONS
! -----

! - Find minimum value in the array (for auto scaling of plots)

FUNCTION minY(inarray() as DOUBLE) AS DOUBLE

 a=inarray(1,2)

 for i=1 to DPi

 if a>inarray(i,2) then

 a=inarray(i,2)

 end if

 next

 minY=a

end function

! - Find maximum value in the array (for auto scaling of plots)

function maxY(inarray() as DOUBLE) AS DOUBLE

 a=inarray(1,2)

 for i=1 to DPi

 if a<inarray(i,2) then

 a=inarray(i,2)

 end if

 next

 maxY=a

end function

! - END OF THE CODE -