

```

! -----
! GRAPHICAL USER INTERFACE (GUI) for CV and SWV measurements
! This is a test program for UWED to make CV and SWV recordings
! and enter parameters on graphical user interface
! Written by Julia Redston, 2016. Whitesides Group, Harvard.
! -----

! -----
! VARIABLES AND OBJECTS
! -----

! -- Labels and text boxes for inputs
dim textBox1 AS TextField , textBox2 AS TextField, textBox3 AS TextField
dim textBox4 AS TextField, textBox5 AS TextField, textBox6 AS TextField
dim label1 AS Label, label2 AS Label, label3 AS Label, label4 AS Label
dim label5 AS Label, label6 AS Label
dim label8 AS Label
dim hzLabel AS Label
dim fileName as TextField
dim startButton AS Button
dim stopButton AS Button
dim settingsButton AS Button
dim firstUni AS INTEGER
firstUni = 0
dim textNote as Label
dim textCycles AS TextField
dim labelCycles AS Label
dim segControl as SegmentedControl
mdone=0
dim dirControl as SegmentedControl
dim cycleControl as SegmentedControl
dim pulseAmp AS DOUBLE
DIM labelPulse AS Label
DIM pulseText AS TextField
DIM startLabel AS Label
DIM startText AS TextField
dim startVol as DOUBLE
Dim stringFile as STRING

! -- Program parameters --
DIM upperVoltage AS Double ! [V] This is upper voltage value
DIM lowerVoltage AS Double ! [V] This is lower voltage value
numberOfCycles = 2! This is number of full cycles
DIM scanRate AS Double ! [mV/s]
DIM timeStep AS Double ! [ms] timestep between measurements divided by 2.

DIM initW AS LONG
initW=22500 ! This is now constant and not calculated
DIM initR AS LONG
initR=22606 ! This should start as close to zero as possible
Rslope=1.5040E4 ! This is number of digits per volt when Ref is changed

DIM numElec !number of electrons

DIM StepHeightAbs AS INTEGER !absolute value of step height
DIM stepR(50) AS INTEGER
DIM nrSteps(50) AS INTEGER
DIM nrDataPoints AS LONG

```

```

DIM nrScans AS INTEGER
dim pulseAmpD AS INTEGER

! Resulting data storage
DIM Vw(10000) AS DOUBLE ! Potential in volts
DIM IuA(10000) AS DOUBLE ! Electrode current in uA
DIM dWe(10000) AS LONG ! Working electrode as set
DIM dRe(10000) AS LONG ! Reference electrode as set
DIM dADC(10000) AS LONG ! ADC value
DIM Ce(10000) AS INTEGER ! Counter electrode potential (as character)
DIM Ni(10000) AS INTEGER ! Index of point in scan
DIM DPi AS INTEGER ! Data point counter
DIM SCi AS INTEGER ! Scan counter
DIM ScanRunning AS INTEGER ! if 0 does not run, if 1 running
ScanRunning=0
DIM lastIndex AS INTEGER
lastIndex=0

DIM TS(10000) AS LONG

! - Graphs visualization
DIM p1 AS Plot ! This is the result plot
DIM pCV AS PlotPoint ! This is one dataset
DIM myCV(10000,2) !holds all cv points
DIM myCV2(10000,2) ! holds all square wave points
DIM graphArray(2000,2) !holds the most recent 2000 points,
! this is what is displayed on graphed

! -----
! SETUP GUI AND APPLICATION
! -----

p1=Graphics.newPlot ! create new plot
pCV=p1.newPlot(graphArray) !set to the Graphing Array
!p1.setTitle("Electrochemical Measurements")
p1.setXaxisLabel("Ref - Work")
p1.setYaxisLabel("Current")

! this allows you to zoom and move around graph, but not click points
p1.setAllowedGestures($F0FF)

p1.setRect( Graphics.width/2 ,Graphics.height/2 - 250 , _
            Graphics.width/2 + 200 ,Graphics.width/2 + 200)
p1.setXaxisLabel("Potential Working vs Ref (V)")
p1.setYaxisLabel("Current working electrode (µA)")
p1.setTitle("Electrochemical Measurements with RFUWED")
textNote = Graphics.newLabel(10, Graphics.height/30,160)

! This forces the app to be locked in the horizontal orientation
! Application shall be launched horizontal
System.setAllowedOrientations(2)

! - initialize and properties of all the buttons and text boxes -

startButton = Graphics.newButton(10,Graphics.height/15)
startButton.setTitle("Start")
startButton.setBackgroundColor(.85, .85, .85)

```

```

stopButton = Graphics.newButton(120,Graphics.height/15)
stopButton.setTitle("Stop")
stopButton.setBackgroundColor(.85, .85, .85)

settingsButton = Graphics.newButton(230,Graphics.height/15)
settingsButton.setTitle("Save Settings")
settingsButton.setBackgroundColor(.85, .85, .85)

segControl = Graphics.newSegmentedControl(10, Graphics.height/15 *2, 200, 30)
segControl.insertSegment("CV",1,0)
segControl.insertSegment("Square Wave",2,0)
segControl.setSelected(1)

label2 = Graphics.newLabel(10, Graphics.height/15 *3,160)
label2.setText("Voltage Upper(V)")

textBox2 = Graphics.newTextField( 150, Graphics.height/15 *3)
textBox2.setBackgroundColor(.85,.85,.85)

label3 = Graphics.newLabel(10, Graphics.height/15 *4,160)
label3.setText("Voltage Lower(V)")

textBox3 = Graphics.newTextField( 150, Graphics.height/15 *4)
textBox3.setBackgroundColor(.85,.85,.85)

label5 = Graphics.newLabel(10, Graphics.height/15 *5,160)
label5.setText("Scan Rate(mV/s)")

textBox5 = Graphics.newTextField( 150, Graphics.height/15 *5)
textBox5.setBackgroundColor(.85,.85,.85)

label6 = Graphics.newLabel(10, Graphics.height/15 *6,160)
label6.setText("Delta Time(ms)")

textBox6 = Graphics.newTextField( 150, Graphics.height/15 *6)
textBox6.setBackgroundColor(.85,.85,.85)

hzLabel = Graphics.newLabel(250, Graphics.height/15 *6,120)
hzLabel.setText("Hz ")

labelCycles = Graphics.newLabel(10, Graphics.height/15 *7,160)
labelCycles.setText("Num Scans")

textCycles = Graphics.newTextField( 150, Graphics.height/15 *7)
textCycles.setBackgroundColor(.85,.85,.85)

label8 = Graphics.newLabel(10, Graphics.height/15 *8,160)
label8.setText("File Name")

fileName = Graphics.newTextField( 150, Graphics.height/15 *8, 180)
fileName.setBackgroundColor(.85,.85,.85)

labelPulse = Graphics.newLabel(10, Graphics.height/15 *9,160)
labelPulse.setText("Pulse Amp(mV)")

pulseText = Graphics.newTextField( 150, Graphics.height/15 *9)
pulseText.setBackgroundColor(.85,.85,.85)
labelPulse.setHidden(1)

```

```

pulseText.setHidden(1)

startLabel = Graphics.newLabel(10, Graphics.height/15 *9,160)
startLabel.setText("Start Value(V)")

startText = Graphics.newTextField( 150, Graphics.height/15 *9)
startText.setBackgroundColor(.85,.85,.85)

dirControl = Graphics.newSegmentedControl(10, Graphics.height/15 *10, 200, 30)
dirControl.insertSegment("Negative ->",2,0)
dirControl.insertSegment("Positive <-",1,0)
dirControl.setSelected(1)

dirControl.setHidden(1)

cycleControl = Graphics.newSegmentedControl(10, Graphics.height/15 *11, 200, 30)
cycleControl.insertSegment("Unidirectional",1,0)
cycleControl.insertSegment("Bidirectional",2,0)
cycleControl.setSelected(1)

cycleControl.setHidden(1)

System.ShowGraphics(1) ! SHOW GRAPHICS

! creating the file settings
! if the settings file exists load the saved settings
if( EXISTS("UWEDSettings.txt")) then
  OPEN "UWEDSettings.txt" FOR INPUT AS #2
  LINE INPUT #2, a$
  if a$ = "Square Wave" then
    segControl.setSelected(2)
    LINE INPUT #2, a$
    textBox2.setText(a$)
    LINE INPUT #2, a$
    textBox3.setText(a$)
    LINE INPUT #2, a$
    textBox5.setText(a$)
    LINE INPUT #2, a$
    textBox6.setText(a$)
    IF NOT (textBox6.getText = "") then
      hzLabel.setText(STR(250/VAL(textBox6.getText))&" Hz")
    end if
    LINE INPUT #2, a$
    textCycles.setText(a$)
    LINE INPUT #2, a$

    pulseText.setText(a$)
    labelPulse.setHidden(0)
    pulseText.setHidden(0)
    dirControl.setHidden(0)
    cycleControl.setHidden(0)
    startLabel.setHidden(1)
    startText.setHidden(1)
  else
    segControl.setSelected(1)
    LINE INPUT #2, a$
    textBox2.setText(a$)
    LINE INPUT #2, a$

```

```

    textBox3.setText(a$)
    LINE INPUT #2, a$
    textBox5.setText(a$)
    LINE INPUT #2, a$
    textBox6.setText(a$)
    IF NOT (textBox6.getText = "") then
        hzLabel.setText(STR(250/VAL(textBox6.getText))&" Hz")
    end if
    LINE INPUT #2, a$
    textCycles.setText(a$)
    LINE INPUT #2, a$
    startText.setText(a$)

    labelPulse.setHidden(1)
    pulseText.setHidden(1)
    dirControl.setHidden(1)
    cycleControl.setHidden(1)
    startLabel.setHidden(0)
    startText.setHidden(0)
end if
CLOSE #2
end if

```

```

! -- BLE --
!starts scanning for ble
DIM bm AS BLEPeripheral
ble.startble
DIM uuid(0) AS STRING
DIM services(1) AS STRING
DIM bleWriteChar AS BLECharacteristic
ble.startscan(uuid)

```

END

```

! -----
! FUNCTIONS
! -----

```

```

!Initialize UWED
! Gets the UWED ready to start, sends constants
SUB InitializeUWED()
    DPi=0 ! No data points yet
    SCi=0 ! Scan counter
    ! Make graph empty
    FOR i=1 TO nrDataPoints
        myCV(i,1)=0
        myCV(i,2)=0
        myCV2(i,1)=0
        myCV2(i,2)=0
    NEXT i
    ! Start setting up
    sendToUWED("A("&STR(initR - Rslope*startVol)&")")
    sendToUWED("E("&STR(timeStep)&")")
    sendToUWED("B("&STR(initW)&")")
    sendToUWED("D(3)")
    sendToUWED("C(1)")
END SUB

```

```

! Initialize a scan
! called by either touchupinside or procdata
! performs 1 scan if scans havent completed
! saves data to file if scans have been completed
SUB InitializeScan()
    SCi=SCi+1 ! Increase scan counter
    system.wait(.1) ! waiting 1/10 a second here prevents skipping data points
    print "on scan ";SCi
    textNote.setText("on scan "&STR(SCi))
    IF SCi>nrScans THEN ! If increase scan counter is larger, DONE
        !save to file
        print "SCi larger than nrScans "
        print SCi; " > "; nrScans
        textNote.setText("writing to file")
        print "writing to file"
        sendToUWED("C(0)") ! Put UWED to potentiometric (off state)
        ! Save data
        if NOT( fileName.getText = "" ) then ! if a file name has been given
            stringFile = fileName.getText & "-raw.txt"
            dim stringFile2 AS String
            stringFile2 = fileName.getText & "-proc.txt"
            OPEN stringFile FOR OUTPUT AS #1 !create raw text file
            OPEN stringFile2 FOR OUTPUT AS #3 !create processed text file
            !first deal with processed data, saves the actual data points graphed
            !aka the data in myCV and myCV2
            IF segControl.title(segControl.selected) = "Square Wave" then
                PRINT #1, "Square Wave "
                PRINT #3, "Square Wave "
                for p=1 to (DPi/2)
                    PRINT #3, STR(myCV2(p,1))&", "STR(myCV2(p,2))!
                next p
            else
                PRINT #1, "CV "
                PRINT #3, "CV "
                print "dpi " &STR(DPi)
                for p=1 to DPi-1
                    PRINT #3, STR(myCV(p,1))&", "STR(myCV(p,2))!
                next p
            end if
            ! raw data
            ! saves reference electrode (dRe), working electrode (DWe), ADC value (dADC)
            ! counter electrode potential (Ce) and index of point (Ni)
            FOR i=1 TO DPi
                PRINT #1, STR(dRe(i))&" "&STR(DWe(i))&" "&STR(dADC(i))&_
                " "&STR(Ce(i))&" "&STR(Ni(i))
            NEXT i
            CLOSE #1
            CLOSE #3
            print "written to file"
        end if
        ! Send email now
        ! if you dont want to send an email comment out all of this
        ! sends from the built in "mail" app and whatever email is saved there
        if NOT( fileName.getText = "" ) then
            ! only sends email with file if there is a file
            DIM em AS Email
            DIM em_mess AS STRING !this is the message of the data
            DIM em_sub AS STRING !this is the subject of the data

```

```

    DIM em_add AS STRING !this is who the email is sent to
    ! IF DATA IS ALWAYS BEING SENT TO SAME ADDRESS
    ! YOU SHOULD EDIT THIS LINE TO ADDRESS BEING SENT TO
    ! example em_add = "johnsmith@gmail.com"
    em_add=""
    em_mess="This is data from UWED measurement"
    em_sub= "UWED Data"
    em=System.newEmail(em_add, em_sub, em_mess)
    ! attaches the processed and raw data files
    em.addAttachment(stringFile, "text/plain")
    em.addAttachment(stringFile2, "text/plain")
    em.send
end if
print "done"
textNote.setText("DONE")
ScanRunning=0
ELSE ! this is not the last scan, send data to UWED
    ! Increase same amount both steps
    print "step height "; stepR(SCi)
    print "numsteps ";nrSteps(SCi)
    if segControl.title(segControl.selected) = "CV" then
        sendToUWED("G("&STR(stepR(SCi))&")")
        sendToUWED("I("&STR(stepR(SCi))&")")
    else
        if cycleControl.title(cycleControl.selected) = "Unidirectional" then
            sendToUWED("A("&STR(initR - Rslope*startVol)&")")
            firstUni =1
        end if
        sendToUWED("G("&STR(0-2*pulseAmpD)&")")
        sendToUWED("I("&STR(2*stepR(SCi)+2*pulseAmpD)&")")
    end if
    sendToUWED("L("&STR(nrSteps(SCi))&")")
    sendToUWED("M()") ! Run the sequence
    ScanRunning=1
END IF
END SUB

! - Process incoming data
! called by BLECharacteristicinfo
! takes in the integer array sent by the UWED
! converts that array into a data point and graphs it
SUB procData(input_array() AS INTEGER)
    PRINT "Index: "&CHR(input_array(1))&CHR(input_array(2))
    ! Feedback, so we can see continuity
    ! First check that there is no duplication
    IF(input_array(1)=lastIndex) THEN
        GOTO OutSub ! If same as last time do nothing
    END IF
    lastIndex=input_array(1) ! Otherwise put it same
    ! Now look if new data came M character on second position
    IF(input_array(2)=ASC("M"))THEN ! New data came
        IF ScanRunning=1 THEN ! Only if scan is running there is something to come
            DPi=DPi+1
            dRe(DPi)=ConvertToNumberBin(input_array,3,4)
            dADC(DPi)=ConvertToNumberBin(input_array,5,6)
            DPi=DPi+1
            dRe(DPi)=ConvertToNumberBin(input_array,7,8)
            dADC(DPi)=ConvertToNumberBin(input_array,9,10)

```

```

Ts(DPi)=ConvertToNumberBin(input_array,11,14)
Ts(DPi-1)=Ts(DPi)
dWe(DPi)=ConvertToNumberBin(input_array,15,16)
dWe(DPi-1)=dWe(DPi)
Ni(DPi)=ConvertToNumberBin(input_array,17,18)
Ni(DPi-1)=Ni(DPi)
! Do the math for visualization
! ODD
VoltW=dRe(DPi-1)*(-6.64898E-5)+1.50310 ! Working electrode pot vs ref
Curr=0.00608*(dADC(DPi-1)-30694)-0.04 ! Current in uA
myCV(DPi-1,1)=VoltW
myCV(DPi-1,2)=Curr
! EVEN
VoltW=dRe(DPi)*(-6.64898E-5)+1.50310 ! Working electrode pot vs ref
Curr=0.00608*(dADC(DPi)-30694)-0.04 ! Current in uA
myCV(DPi,1)=VoltW
myCV(DPi,2)=Curr
if segControl.title(segControl.selected) = "CV" then ! CV calculations
    graphArray((DPi MOD 2000 )+ 1 , 1) = myCV(DPi,1)
    graphArray((DPi MOD 2000 )+ 1 , 2) = myCV(DPi,2)
    graphArray(((DPi -1) MOD 2000 )+ 1 , 1) = myCV(DPi-1,1)
    graphArray(((DPi -1)MOD 2000 )+ 1 , 2) = myCV(DPi-1,2)
    ! Show graph
    if( DPi MOD 3 = 0 ) then ! only update graph every 3 points
        pCV.setPoints(graphArray)
    end if
else ! Square wave calculations
    ! if its the first data point in a unidirectional sweep
    ! then set points to 0
    ! prevents noise
    if firstUni = 1 then
        myCV2(DPi/2,1) = (myCV(DPi,1)+myCV(DPi-1,1))/2
        myCV2(DPi/2,2) = 0
        firstUni = 0
    else ! if its not the first point of unidirectional
        ! then perform calculations for square wave
        myCV2(DPi/2,1)=(myCV(DPi,1)+myCV(DPi-1,1))/2 ! average two
        myCV2(DPi/2,2)=myCV(DPi,2)-myCV(DPi-1,2) ! different between two
        graphArray((DPi/2 MOD 2000 )+ 1 , 1) = myCV2(DPi/2,1)
        ! put into the graphing array
        graphArray((DPi/2 MOD 2000 )+ 1 , 2) = myCV2(DPi/2,2)
        print myCV2(DPi/2,1);", "; myCV2(DPi/2,2)
        ! print "dpi "; dpi
        print "dpi mod ";( (DPi/2 MOD 2000) +1)
    end if
    ! Show graph
    if( DPi MOD 3 = 0) then ! only update graph every 3 points
        pCV.setPoints(graphArray)
    end if
end if
END IF
END IF
! If end of the scan start new
IF (input_array(2)=ASC("m")) THEN ! End of scan
    InitializeScan ! Start next scan
END IF
! In other cases do currently nothing
OutSub:

```


END SUB

! - Convert to Number

! called by procdata converts integer array into characters then those chars into a long

FUNCTION ConvertToNumber(input_array() AS INTEGER, _

begi AS INTEGER, _

endi AS INTEGER) AS LONG

 DIM nets AS STRING

 nets=""

 FOR i=begi TO endi

 nets=nets & chr(input_array(i))

 NEXT i

 ConvertToNumber=VAL(nets)

END FUNCTION

! -----

! EVENTS

! -----

! Called, when button is clicked

! Takes in the button that was clicked and time

SUB touchUpInside(ctrl AS Button, when AS DOUBLE)

 mdone=0

 DIM writeValue(1) AS INTEGER

 if NOT (bm = NULL) AND NOT (bleWriteChar = NULL) then ! if ble is connected

 if ctrl = stopButton then ! if stop button send "STOP" command to UWED

 sendToUWED("N()")

 else if ctrl = settingsButton then

 ! if its the settings button, save current settings

 ! save stuff to file

 print "writing settings file"

 OPEN "UWEDSettings.txt" FOR OUTPUT AS #2

 if segControl.title(segControl.selected) = "Square Wave" then

 PRINT #2, "Square Wave"

 PRINT #2, textBox2.getText ! upper voltage

 PRINT #2, textBox3.getText ! lower voltage

 PRINT #2, textBox5.getText !scan rate

 PRINT #2, textBox6.getText !time Step

 PRINT #2, textCycles.getText ! num cycles

 PRINT #2, pulseText.getText ! pulse amp

 else

 PRINT #2, "CV"

 PRINT #2, textBox2.getText ! upper voltage

 PRINT #2, textBox3.getText ! lower voltage

 PRINT #2, textBox5.getText !scan rate

 PRINT #2, textBox6.getText !time Step

 PRINT #2, textCycles.getText ! num cycles

 PRINT #2, startText.getText ! start vol

 CLOSE #2

 end if

 else if ctrl = startButton then

 ! if start button, make calculations start first scan

 DPi =0

 ! clear graphing array

 for p= 1 to 2000

 graphArray(p,1) =0

 graphArray(p,2)=0

 NEXT P

```

upperVoltage= VAL(textBox2.getText) ! [V] This is upper voltage value
lowerVoltage= VAL(textBox3.getText) ! [V] This is lower voltage value
scanRate=VAL(textBox5.getText )      ! [mV/s] change in voltage per second
timeStep=VAL(textBox6.getText)
    ! [ms] time between each step (2x for cv 4x for square)
nrScans = VAL(textCycles.getText)    ! number of scans
numberOfCycles = (nrScans -1 )/2    ! number of cycles
print "numscans total ";nrScans
nrDataPoints=2*numberOfCycles*(nrSteps2+1)
p1.setView(lowerVoltage, -50, upperVoltage, 50, 0)

```

```

IF segControl.title(segControl.selected) = "Square Wave" then ! if square wave
    ! calculate number of steps
    nrSteps1=1000000*ABS(upperVoltage - lowerVoltage)/(scanRate*timeStep*2)
    ! calculate stepheight
    stepHeightAbs=ABS((upperVoltage-lowerVoltage)*Rslope)/nrSteps1
    ! calculate number data points
    nrDataPoints=nrSteps1+2

```

```

    IF(dirControl.title(dirControl.selected) = "Negative ->")THEN ! negative
        stepR(1)=(-1)*stepHeightAbs !starts in negative direction
        startVol = lowerVoltage !starts at lower voltage
    ELSE ! positive
        stepR(1)=stepHeightAbs ! starts positive direction
        startVol = upperVoltage !starts at the upper voltage
    END IF
    nrSteps(1)=nrSteps1
    nrDataPoints=nrSteps1+2
    pulseAmp = VAL(pulseText.getText)
    pulseAmpD = (pulseAmp/1000)*Rslope !Pulse amplitude
    DIM myCV2(nrDataPoints)

```

```

    dim directionalVec as INTEGER
        ! controls if switch directions every scan or not
    directionalVec = -1
    if cycleControl.title(cycleControl.selected) = "Unidirectional" then
        directionalVec = 1
    end if
    print "dirVec ";directionalVec
    for i = 1 TO nrScans
        nrSteps(i) = nrSteps1
        if i MOD 2 = 0 then
            stepR(i) = (directionalVec * stepR(1))
            print "stepR at ";i;" ";stepR(i)
        else
            stepR(i) = stepR(1)
        end if
    NEXT i
else ! CV calculations
    DIM myCV(nrDataPoints)
    startVol = 0
    if(NOT (startText.getText = "")) then !if start text is not empty
        startVol = VAL(startText.getText)
    end if
    ! calculate num steps from 0 to upper voltage
    nrSteps1=1000000*(upperVoltage)/(scanRate*timeStep*2)
    dim nrSteps1Start as double
    ! calculate num steps from start voltage to upper voltage

```



```

IF peripheral.bleName="UWED" THEN ! Right device found
    bm=peripheral
    ble.connect(bm)
    ble.stopScan
    print "Device is UWED"
END IF
END SUB

! - When information about the device is fetched
SUB BLEPeripheralInfo(  time AS DOUBLE, _
                        peripheral AS BLEPeripheral, _
                        kind AS INTEGER, _
                        message AS STRING, _
                        err AS LONG)

IF kind=1 THEN ! Connection completed
    peripheral.discoverServices(uuid)
ELSE IF kind=2 OR kind=3 THEN ! Connection lost
    ble.connect(bm)
ELSE IF kind=4 THEN ! Service dound
    DIM avServ(1) AS bleservice
    avServ=peripheral.services
    FOR a=1 TO UBOUND(services,1)
        FOR a=1 TO UBOUND(services,1)
            IF avServ(a).uuid="2220" THEN ! If right service id has been found
                peripheral.discoverCharacteristics(uuid, avServ(a))
                print "Service discovered"
            END IF
        NEXT
    NEXT
END IF
END SUB

! - When information about service is obtained
SUB BLEServiceInfo(  time AS DOUBLE, _
                    peripheral AS BLEPeripheral, _
                    service AS BLEService, _
                    kind AS INTEGER, _
                    message AS STRING, _
                    err AS LONG)

IF kind=1 THEN
    DIM chx(1) AS blecharacteristic
    chx=service.characteristics
    FOR i=1 TO UBOUND(chx,1)
        IF service.uuid="2220" THEN
            IF chx(i).uuid="2221" THEN
                peripheral.setNotify(chx(i),1) ! Notify if changes
                print "Read characteristic found"
            ELSE IF chx(i).uuid="2222" THEN
                bleWriteChar=chx(i)
                print "Write characteristic found"
            END IF
            textNote.setText("BLE CONNECTED!")
        END IF
    NEXT
END IF
END SUB

! - If BLE Characteristic info comes

```

```

! when UWED sends back data
SUB BLECharacteristicInfo(      time AS DOUBLE, _
                                peripheral AS BLEPeripheral, _
                                characteristic AS BLECharacteristic, _
                                kind AS INTEGER, _
                                message AS STRING, _
                                err AS LONG)

    IF kind=2 THEN
        IF characteristic.uuid="2221" THEN
            DIM value(1) AS INTEGER
            value=characteristic.value
            procData(value)
        END IF
    END IF
END SUB

! -----
! HELPING FUNCTIONS
! -----

! currently not used
! could be used to autoscale the graph, but it makes the program run much slower
FUNCTION minY(inarray() AS DOUBLE) AS DOUBLE
    a=inarray(1,2)
    for i=1 to DPi
        if a>inarray(i,2) then
            a=inarray(i,2)
        end if
    next
    minY=a
end function

! currently not used
! could be used to autoscale the graph, but it makes the program run much slower
function maxY(inarray() AS DOUBLE) AS DOUBLE
    a=inarray(1,2)
    for i=1 to DPi
        if a<inarray(i,2) then
            a=inarray(i,2)
        end if
    next
    maxY=a
end function

! - Convert to Number from String source - currently not used
FUNCTION ConvertToNumberStr(input_array() AS INTEGER, _
                            begi AS INTEGER, _
                            endi AS INTEGER) AS LONG

    DIM nets AS STRING
    nets=""
    FOR i=begi TO endi
        nets=nets & chr(input_array(i))
    NEXT i
    ConvertToNumberStr=VAL(nets)
END FUNCTION

! - Convert to Number from Binary source - used in data receiving
FUNCTION ConvertToNumberBin(input_array() AS INTEGER, _

```

```

                                begi AS INTEGER, _
                                endi AS INTEGER) AS LONG

DIM sum AS LONG
sum=input_array(begi)
FOR i=begi+1 TO endi
    sum=sum*256
    sum=sum+input_array(i)
NEXT i
ConvertToNumberBin=sum
END FUNCTION

! - when the segcontrol's value is changed (when a tab is clicked)
SUB valueChanged ( ctrl AS Control, _
                    time AS DOUBLE)
IF ctrl = segControl then ! change segControl value
    ! if set to CV, make appropriate textboxes etc visible
    IF segControl.title(segControl.selected) = "CV" then
        labelPulse.setHidden(1)
        pulseText.setHidden(1)
        dirControl.setHidden(1)
        cycleControl.setHidden(1)
        startLabel.setHidden(0)
        startText.setHidden(0)
    Else
        ! if set to Square Wave, make appropriate textboxes etc visible
        labelPulse.setHidden(0)
        pulseText.setHidden(0)
        dirControl.setHidden(0)
        cycleControl.setHidden(0)
        startLabel.setHidden(1)
        startText.setHidden(1)
    END IF
ELSE IF ctrl = textBox6 then
    hzLabel.setText((250/VAL(textBox6.getText))&" Hz")
END IF
END SUB

! - END OF THE CODE -

```